

Real-Time Instrument Control of the Orbitrap Tribrid Mass Spectrometer

Derek J. Bailey, Florian Grosse-Coosmann, Manish Doshi, Qingyu Song, Jesse D. Canterbury, Qiming Wan, and Michael W. Senko, Thermo Fisher Scientific, 355 River Oaks Parkway, San Jose, CA 95134

OVERVIEW

Presented here is an instrument application programming interface (IAP) for real-time control of Thermo Scientific™ Orbitrap™ Tribrid™ mass spectrometers (MS). The IAP provides a programmatic and interactive way to receive and send data between the MS and host computer system.

INTRODUCTION

In the past half-decade a growing emphasis has been placed on the ability to incorporate real-time data analysis and programmatic control during the mass spectrometer acquisition. To provide advanced, high-performance access for third-parties, we previously released an IAP for the bench-top Orbitrap mass spectrometers (Exactive-series IAP) (ES-IAP). We have now enabled the Orbitrap Tribrid family of instruments to use a similar, expanded programming interface to extend real-time capabilities to researchers on additional platforms.

METHODS

The Tribrid IAP is written for the Microsoft.NET Framework (Version 4.5.1) and is fully integrated with the Tune instrument control software (Version 2.2). It uses an event-driven architecture where users subscribe to instrument generated events throughout the acquisition (e.g., spectra, messages, device readbacks, etc.). Users are able to control the instrument in real-time by modifying instrument parameters or submitting custom scan definitions. The IAP is efficient and does not impact the performance of the instrument. The Tribrid IAP is a direct extension of the previous bench-top Orbitrap interface, but various improvements and refinements break binary compatibility. All example code and applications are written in C#.

RESULTS

Figure 1. Compatible Mass Spectrometers




Figure 2. IAP Architecture




Figure 3. IAP Instance Creation Example

```
IFusionInstrumentAccessContainer iAccessContainer  
= Factory<IFusionInstrumentAccessContainer>.Create();
```

Creating an instance of the Tribrid IAP is straightforward using the factory method pattern provided by the *Thermo.TNG.Factory.dll* assembly (Figure 3). There is no longer a need of using COM or manually searching the registry to create an instance from an explicitly loaded assembly. All other access to the IAP can be obtained through the factory-generated *IFusionInstrumentAccessContainer* instance.

SPECTRUM ACCESS

The spectrum-related interfaces are broken into 5 separate parts as outlined in Figure 4.

- IMassIntensity**
Represents a m/z-intensity pair
- INoiseNode**
Represents a noise baseline of a spectrum. It derives from *IMassIntensity*
- ICentroid**
Represents the instrument-centroided peak and associated data. If profile data is present, an array of *IMassIntensity* will represent those points. It derives from *IMassIntensity*.
- ISpectrum**
Represents a mass spectrum of *ICentroid* peaks and its associated polygonal noise-baseline. The mass analyzer that collected the data is indicated by the detector name (e.g., Ion Trap or Orbitrap).
- InformationSourceAccess**
A dictionary-like container to store various metadata fields (e.g., Scan Trailer, Status Log, Tune Data, etc.) that are readily present in the Thermo raw data file. Provides access to both the string representation as well as the raw-values (.NET object) for each item.


Figure 5. Example Spectrum Access

```
private void ExampleSpectrum(ISpectrum spectrum)  
{  
    // Check for emptiness  
    if (!spectrum.CentroidCount.HasValue) return;  
  
    // Loop over the spectrum centroids  
    foreach(ICentroid centroid in spectrum.Centroids)  
    {  
        // Print some information about each peak to console  
        Console.WriteLine("MZ: {0} Intensity: {1} Charge: {2}",  
            centroid.Mz,  
            centroid.Intensity,  
            centroid.Charge.HasValue ? centroid.Charge.Value.ToString() : "?");  
    }  
}
```

EVENTS


The IAP uses an event-driven model for informing the user application of instrument-related events. When the instrument has finish collecting a spectrum, it publishes the event and listening clients can consume it (Figure 6). The same spectrum event the IAP uses is what the Tune application and the raw data file internally subscribe to; this ensures the data displayed in Tune and stored in the raw data file is 100% identical to the one the IAP user would receive.

Figure 6. Spectrum Receiving Interfaces



The spectrum receiving event for the IAP is controlled through the *IMsScanContainer* interface (Figure 7.).

Figure 7. Spectrum Receiving Interfaces



Instrument Control

The IAP allows users to submit scans, at will, for the instrument to perform. This capability is available all the time—when running a method or controlled by Tune. There are two types of user-submitted scans: **Repeating Scans** and **Custom Scans**.

Figure 8. Instrument Scan Priority Stack





Figure 9. Example IAP scans During Method Acquisition




A top speed (3s) data-dependent method was acquired while an IAP application sent both repeating and custom scans. The method logic was suspended until the IAP scan events were completed and then resumed where it left off. IAP scans inserted during method acquisition are placed in analysis pipeline so that they can be parallelized for optimal performance. This insertion may lead to a small delay between when a scan is submitted and its acquisition.

INSTRUMENT CONTROL

The IAP offers a convenient way for specifying the parameters of a scan. After creating a default repeating or custom scan through an IAP method call, simply specify the parameters in a *Dictionary<string,string>* Values property of the scan (Figure 10). The allowable keys and values ranges for the scans are provided by another IAP property that is not shown here.

Figure 10. Example Energy Ramp of a Precursor Recorded from Tune




While acquiring spectra from Tune, the above code was executed from an IAP application. Here 5 custom scans were created, targeting a peak at 524.3 m/z for MS/MS analysis. This peak was fragmented between 20 and 40 normalized collision energy (NCE) in increments of 5. These custom scans were immediately executed and saved in the acquiring raw data file. After these custom scans were acquired, the Tune-defined scan continued for the rest of the acquisition.

EXAMPLE IAP APPLICATIONS

Spectral Matching in Real Time



The Tribrid Method Editor offers complex experimental design and real-time filtering of mass spectrum, but is limited to the built-in filters and execution workflow. The IAP offers a way for researchers to extend the capabilities by implementing their own decision-making logic and control. Here we demonstrate real-time spectral matching to decide whether to take an MS3 scan.

Figure 11. Spectral Library Creation



A standard MS/MS experiment was constructed in the Method Editor application and acquired on the Fusion mass spectrometer. Some of the resulting MS/MS spectra were extracted from the raw data file and imported in a local SQLite database. The spectra were binned into single Dalton bins and intensities were normalized from 0 to 100.



Figure 12. Spectral Matching Application



Real Time Compound Collision Energy Optimization

In typical data-dependent (DD) experiments, the MS/MS collision energy (CE) is set to a static value that performs well over a broad range of precursors. In order to better optimize the collision energy on a precursor-by-precursor basis, we wrote a simple IAP application that performs a series of rapid ion trap scans over a CE range. The application then determined the optimal CE and subsequently performed a high resolution Orbitrap scan of the same precursor at the optimal CE (Figure 11.).

Figure 11. Compound Collision Energy Optimization



While acquiring spectra from Tune, both 262.6 and 524.3 m/z were targeted for MS/MS analysis. A series of rapid ion trap (IT) scans were performed on each precursor at increasing NCEs. Once the TIC of the of the IT survey scan was less than 1% of the TIC at 0 NCE, the optimization was stopped and a high resolution Orbitrap MS/MS analysis was queued at the optimal NCE. Both of these compounds were optimized in tandem and arrived at different optimal NCEs. The 262 precursor broke apart quicker than the 524 precursor.

CONCLUSIONS

The Instrument Application Programming Interface (IAP) for the Orbitrap Tribrid Mass Spectrometer family provides advanced control and acquisition to enable custom experimentation and method development for mass spectrometry. This increased flexibility the IAP provides users has the potential to spawn creative new experimental design and experiments not possible with the standard software.

The IAP offers a lot of functionality not mentioned in this document. For complete listing of the capabilities, as well as all the example applications shown, please visit our Thermo Fisher Scientific Life Science Mass Spectrometry GitHub page: <https://github.com/thermofisherlscms/iapi>

ACKNOWLEDGEMENTS

We thank Phil Remes for his suggestions and technical support implementing the IAP on the Tribrid system. We also thank Baram Barange for contributing to the example programs.

TRADEMARKS/LICENSING

Access to the API is not part of the standard software delivered with the instrument and requires a special license agreement. Contact derek.bailey@thermofisher.com for details.

© 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.